









**BI** **NA** **RE** **E**  
make **visible**

TISSUE CLEARING

PROTOCOL



NAME OF PROTOCOL

**The optical clearing protocol for Spheroid**


Cat.No. HRSC-101

CONDITION OF SAMPLE : Spheroid  
≤ Length 1 mm x Width 1 mm x Height 1 mm

CODE OF PROTOCOL : C1001

REVISION OF PROTOCOL : 1.1.9 (2021.01.02)

[A] - Preparation | Planning you test



[Spheroid Clearing & Imaging within 2 days]

When we designed the protocol, we considered not only the effectiveness of the clearing but also the working time of the researchers.

Enjoy the tissue clearing!

[B] - Preparation | Taking the solutions

B-1. All of the solutions should be stored at 4°C.

B-2. Check Spheroid Clearing Solution and Mounting & Storage Solution for crystallization or precipitation before each use.  
Redissolve any precipitation by warming the solution at 37°C for 1-2 h and then use.

B-3. Do not use the individual solutions from the other kit. Even if the names of solutions are the same.

The component compositions are not the same. Each solution has a unique component composition depending on the purpose of the kit.

① Starting Solution

② Spheroid Clearing Solution

+ Mounting & Storage Solution

- The Mounting & Storage Solution (Cat. No. SHMS-060) is not included in Binaree Tissue Clearing Kit for Spheroid (HRSC-001). Only the Starter's kits (HRSC-101) contain the Mounting & Storage Solution.

- The solutions may become crystallized or precipitated. If this occurs, incubate it at 37°C for 1-2 h before use.

[C] - Preparation | Fixing the sample

C-1. Incubate the sample with 4% PFA at 4°C for 15 min.

C-2. Wash the sample with 1 x PBS while shaking at 4°C for 10 min X 3 times.

C-3. Incubate the ② Spheroid Clearing Solution and + Mounting & Storage Solution at 37°C for 1-2 h before use.

## [D] - Protocol | Clearing the fixed sample

- D-1. Incubate the sample with ① 0.5 ml Starting Solution at 4° C until the sample sinks.
- D-2. Incubate the sample with ② 0.5 ml Spheroid Clearing Solution in a shaking at 50 rpm /37°C for 12-24 h.
- D-3. Wash the sample with **distilled water** while shaking at 50 rpm/ 4°C for 10 min X 3 times.  
The sample may become opaque and swell. This does not affect the clearing process --> The sample will be cleared again in Mounting & Storage Solution.  
If the sample not enough clear in D-2, spheroid clearing (D-2) & washing (D-3) should be repeated until cleared.
- D-4. (optional) Add nuclear stain solution (e.g. DAPI, 20-40 µg/ml in **distilled water**) while shaking at 4°C for 1-2 h.
- D-5. Incubate the sample with + 0.3 ml Mounting & Storage Solution in a shaking incubator at 50 rpm/ 37°C for at least 12 h or more.

## [E] - Clearing Tips

- E-1. If the sample contains air bubbles → Centrifuge the sample at 3,000 rpm/24°C for 1 min.
- E-2. If the sample is not entirely cleared → Repeat from step D-2 to step D-3.
- E-3. If the rpm is not specified → Operate the shaking incubator gently.
- E-4. **Never wash the sample with PBS** instead of distilled water at step D-3.
- E-5. It is recommended to use the vial for clearing rather than the chamber slide. Taking image via confocal microscope, use a slide chamber (2 wells or 4 wells) like the image below. Sealing the chamber with label tape reduces drying.




Figure 1. When taking images through confocal microscopy, the image chamber must be seal by label tape.

## [F] - Storage & Imaging Tips

- F-1. Store the cleared sample in +Mounting & Storage Solution at the room temperature (20~25°C).
- F-2. **Take images within 3 days after the clearing** for the best results.
- F-3. Take images on the microscope. We recommend using a Confocal Laser Scanning Microscope (CLSM).
- F-4. +Mounting & Storage Solution is a **solvent-free** material that is safe to use in the Light Sheet Fluorescence Microscope (LSFM).
- F-5. **Refractive Index(RI)** of the +Mounting & Storage Solution is 1.46.
- F-6. Be careful of making bubbles while filling the microscope chamber with the sample and the +Mounting & Storage Solution.  
The bubbles may disturb the imaging.

## [G] - Contact Us | Technical support

Binaree, Inc. (Headquaters)

- 47 Gyeongdaero17-gil Buk-gu, STE#608 IT Convergence Bldg(115)., Daegu, 41566, Republic of Korea.
- Website: binaree.com ○ Email : lab@binaree.com
- Tel : +82-(0)53-939-5012 ○ Fax : +82-(0)53-382-5012