




TISSUE CLEARING

PROTOCOL

HRTC-002




NAME OF PROTOCOL

Binaree Tissue Clearing™ kit (High Resolution)

Cat. No. HRTC-002

CODE OF PROTOCOL: C1001  
REVISION OF PROTOCOL: 1.1.8 (2022.08.17)

## 1. Preparation | Planning your test



When we designed the protocol, it was considered that clearing effectiveness and working time are the most important factors in the protocol process.

Enjoy the tissue clearing!

## 2. Preparation | Taking the solutions

2-1. All the solutions should be stored at 4 °C.

2-2. Check Tissue Clearing Solution A & B and Mounting & Storage Solution for crystallization or precipitation before use.

2-3. Contents of the kits

	Binaree Tissue Clearing kit (HRTC-002)	Mounting & Storage Solution (SHMS-060)
① Starting Solution	30 ml	
② Tissue Clearing Solution A	30 ml	
③ Tissue Clearing Solution B	60 ml	
+ Mounting & Storage Solution		60 ml

- The Mounting & Storage Solution (Cat. No. SHMS-060) is not included in Binaree Tissue Clearing Kit (HRTC-002).
- The solutions may become crystallized or precipitated. If this occurs, incubate it at 37 °C for 1-2 h before use.

## 3. Preparation | Fixing the sample

3-1. The mouse is transcardial perfused with 4% PFA.

3-2. Incubate the sample with 4% PFA at 4 °C for overnight.

3-3. Wash the sample with 1X PBS while shaking at 4 °C for 20 min × 3 times.

**Note.** Incubate the ②③ Tissue Clearing Solution A & B and + Mounting & Storage Solution at 37 °C for 1-2 h before use.

NAME OF PROTOCOL

**Binaree Tissue Clearing™ kit (High Resolution)**

Cat. No. HRTC-002

#### 4. Protocol I Clearing the fixed sample

Please refer to APPENDIX and adjust the clearing step according to the tissue thickness.

4-1. Incubate the sample with ① **10 ml Starting Solution** at 4 °C until the sample sinks.

4-2. Transfer the sample to ② **10 ml Tissue Clearing Solution A** and agitate at 50 rpm /37 °C for 48 h.

4-3. Wash the sample with distilled water while shaking at 50 rpm /4 °C for 30 min × 4 times.

The sample may become opaque and swell. This does not affect the clearing process; the sample will be cleared again in Mounting & Storage Solution.

4-4. Transfer the sample to ③ **10 ml Tissue Clearing Solution B** and agitate at 50 rpm /37 °C for 48 h.

If the tissue does not enough clear in step 4, washing (4-3) & tissue clearing (4-4) should be repeated until cleared.

4-5. Wash the sample with distilled water while shaking at 50 rpm /4 °C for 30 min × 4 times.

The sample may become opaque and swell. This does not affect the clearing process; the sample will be cleared again in Mounting & Storage Solution.

4-6. (optional) Add nuclear stain solution (e.g. DAPI, 20-40 µg/ml in distilled water) while shaking at 4 °C for overnight.

Wash the sample with distilled water while shaking at 50 rpm /4 °C for 30 min × 3 times.

4-7. Transfer the sample to + **20 ml Mounting & Storage Solution** and agitate at 50 rpm /37 °C for 12-24 h.

#### 5. Clearing Tips

5-1. If the sample contains air bubbles → Centrifuge the sample at 3,000 rpm / 24 °C for 1 min.

5-2. If the sample is not entirely cleared → Repeat from step 4-3 to step 4-4.

5-3. If the rpm is not specified → Operate the shaking incubator gently.

5-4. In steps 4-3 and 4-5, the sample may be washed with 0.1X PBS instead of distilled water.

5-5. It is recommended to use the vial for tissue clearing rather than the chamber slide.

Drying causes crystallization of Tissue Clearing Solution A & B and Mounting & Storage Solution.

#### 6. Storage & Imaging Tips

6-1. After imaging, the samples are stored at room temperature in the Mounting & Storage Solution. Never refrigerate.

6-2. Take images within 7 days after the clearing for the best results.

6-3. Take images on the microscope. We recommend using a Light Sheet Fluorescence Microscope (LSFM) or Confocal Laser Scanning Microscope (CLSM). Analyze and visualize the images with a microscopy image analysis software.

6-4. Mounting & Storage Solution is a solvent-free material that is safe to use in the Light Sheet Fluorescence Microscope (LSFM).

6-5. Refractive Index (RI) of the Mounting & Storage Solution is **1.45**.

NAME OF PROTOCOL

Binaree Tissue Clearing™ kit (High Resolution)

Cat. No. HRTC-002

6-6. Be careful of making bubbles while filling the microscope chamber with the sample and the Mounting & Storage Solution. The bubbles may disturb the imaging.

6-7. To take images of tissue with less than 1 mm thickness via confocal microscope, use a slide chamber (2 wells or 4 wells) like the image below. Sealing the chamber with label tape reduces drying. Too much of the Mounting & Storage Solution can cause the sample in the chamber to shake. The optimal volume for 1 mm-thick tissue is 200 µl.




Figure 1. When taking images through confocal microscopy, the image chamber must be seal by label tape.

## 7. Appendix

Read not only the appendix but also the protocol. The protocol describes the method in detail.

### 7-1. Mouse ≤ 8-weeks-old

step	Summary	Temp. (°C)	1 mm thickness	3 mm thickness	≤ 7 mm thickness (ex. half brain, spleen, spinal cord)	≥ 7 mm thickness (ex. whole brain, lung, kidney, heart)	Spheroid Organoid
3-2	4% PFA	4 °C	12 - 24 h	12 - 24 h	12 - 24 h	12 - 24 h	15 min - 1 h
3-3	Wash with 1X PBS	4 °C	20 min × 3 times	20 min × 3 times	20 min × 3 times	20 min × 3 times	20 min × 3 times
4-1	Starting Solution*	4 °C	2 - 12 h	12 - 24 h	24 h	24 - 48 h	2 - 12 h
4-2	Tissue Clearing Solution A	37 °C	4 - 6 h	24 - 48 h	48 h	48 - 72 h	
4-3	Wash with distilled water	4 °C	10 min × 3 times	30 min × 4 times	30 min × 4 times	30 min × 4 times	
4-4	Tissue Clearing Solution B	37 °C		24 - 48 h	24 - 48 h	48 - 72h	4 - 24 h
4-5	Wash with distilled water	4 °C		30 min × 4 times	30 min × 4 times	30 min × 4 times	10 min × 4 times
4-6	(option) DAPI staining	4 °C	1 h	6 - 12 h	12 - 24 h	12 - 24 h	1 h
4-7	Mounting & Storage Solution	37 °C	< 1 day	> 1 day	> 1 day	> 2 days	> 6 h

**Note.** \*Samples were incubated in Starting Solution until the sample sank. However, even if the sample (ex. lung, spinal cord, spheroid, etc.) does not sink in the Starting Solution after 2 days, proceed to the next steps.

**Note.** If the tissue does not enough clear in step 4-4, washing (step 4-3) & tissue clearing (step 4-4) should be repeated until cleared.

**Note.** When using tissue from mice older than 8 weeks, repeat the tissue clearing (step 4-4) & washing (step 4-5) until the tissue is clear.

**Note.** Spheroid and organoid are not incubated with Tissue Clearing Solution A.

NAME OF PROTOCOL

Binaree Tissue Clearing™ kit (High Resolution)

Cat. No. HRTC-002

## 8. Protocol I (Option) Immunostaining for 1 mm-thick tissue

### 8-1. Tissue Preparation

- 8-1-1. Perfuse mouse with 4% PFA using transcardial perfusion.
- 8-1-2. Post-fix sample with 4% PFA overnight (12-16 h) at 4 °C. Do not over-fix the sample.
- 8-1-3. Wash the sample with 1X PBS while shaking at 4 °C for 20 min × 3 times.
- 8-1-4. Cut tissues into 1 mm-thick sections using Mouse Brain Slicer.

### 8-2. Tissue Clearing (1 mm-thick tissue)

- 8-2-1. Incubate the sample with 6 ml Starting Solution while shaking at 4 °C until the sample sinks (at least 24 h or more).
- 8-2-2. Incubate the sample with 3 ml Tissue Clearing Solution A in a shaking incubator at 37 °C for 4-6 h.
- 8-2-3. Wash the sample with distilled water while shaking at 4 °C for 30 min × 4 times.

**Note.** Tissue may become opaque and swell at step 8-2-3.

If the tissue does not enough clear in step 8-2-2, washing (8-2-3) & tissue clearing solution B should be used to incubate until cleared (at least 12 h or more).

### 8-3. Permeabilization and Blocking

You can perform standard immunostaining procedures.

- 8-3-1. Incubate the sample with permeabilization solution in a shaking incubator at 37 °C for 2-3 days.

**Note.** You can perform standard immunostaining procedures. Permeabilization is a key step for successful antibody labeling. Therefore, the permeabilization condition should be optimized depending on the tissue and antibodies.

**Note.** Sample may become opaque at permeabilization step. Permeabilization solution (recommendation):

1X PBS or TBS containing 0.2-0.5% Triton X-100, 10-20% DMSO, Serum.

### 8-4. Antibody labeling

- 8-4-1. Incubate the sample with primary antibody in a shaking incubator at 37 °C for 3-6 days.
- 8-4-2. Wash the sample with 1X PBS or TBS while shaking at 4 °C for 20 min × 6 times.
- 8-4-3. Incubate the sample with secondary antibody in a shaking incubator at 37 °C for 3-6 days.
- 8-4-4. Wash the sample with 1X PBS or TBS while shaking at 4 °C for 20 min × 6 times.
- 8-4-5. (optional) Add nuclear stain solution (e.g. DAPI, 20-40 µg/ml in distilled water) while shaking at 4 °C for 1 h.
- 8-4-6. Wash the sample with 1X PBS or TBS while shaking at 4 °C for 10 min × 3 times.
- 8-4-7. Incubate the sample with 5 ml Mounting & Storage Solution in a shaking incubator at 37 °C for at least 12 h or more.
- 8-4-8. Image on a Confocal Laser Scanning Microscope.

NAME OF PROTOCOL

Binaree Tissue Clearing™ kit (High Resolution)

Cat. No. HRTC-002

**Note.** The specific depth of penetration for each antibody will be dependent upon the tissue type, the dilution, as well as the incubation time.

**Note.** Successful immunostaining is strongly dependent on the primary and secondary antibodies.

The followings are guidelines for successful immunostaining.

- a. Antibody concentration should be optimized: Primary antibody dilution of 1 : 50 to 1 : 200 and secondary antibody dilution of 1 : 100 to 1 : 500 are required.
- b. Usually, a 1 : 100 dilution of the secondary antibody works well.
- c. Antibody incubation time: for 1 mm-thick tissue, 3-6 days are required.
- d. Antibody dilution buffer (recommendation): 1 X PBS or TBS containing 0.2-0.5% Tween 20, 5% DMSO and 5% serum.
- e. In order to reduce the total immunostaining time, fluorescence conjugated 1st antibody can be used.
- f. In order to reduce the incubation time for the secondary antibody, you can also use an immunoglobulin fragment [e.g. Fab and F(ab')<sub>2</sub>].
- g. Before proceeding to a thicker tissue, immunostaining should be optimized starting with a thinner tissue.
- h. To take images of tissue with less than 1 mm thickness via confocal microscope, use a slide chamber (2 wells or 4 wells) like the image below. Sealing the chamber with label tape reduces drying. Too much of the Mounting & Storage Solution can cause the sample in the chamber to shake. The optimal volume for 1 mm-thick tissue is 200 µl.




Figure 2. When taking images through confocal microscopy, the image chamber must be seal by label tape.

- i. During imaging, the temperature should be above 25 °C. At low temperatures, crystallization or precipitation can form.

## 9. Contact Us | Technical support

Binaree, Inc. (Headquarters)

- Address: 47 Gyeongdaero17-gil Buk-gu, STE#608 IT Convergence Bldg. (115), Daegu, 41566, Republic of Korea.
- Website: binaree.com
- Email: lab@binaree.com
- Tel: +82-(0)53-939-5012
- Fax: +82-(0)53-382-5012